



- Integrated safety lid for secure storage
- Stackable design with optimized airflow
- 100% aluminum & stainless steel construction, no plastic
- Heat-resistant silicone-insulated handles
- Six color-coded lids for easy ID
- Compatible with valve, PTFE, and barrier filter cassettes
- Valve and barrier filter designed for near-lifelong use
- PTFE cassette supports up to 5,200 sterilization cycles
- Structurally enhanced walls retain shape through repeated high-temperature exposure



Valve System



Barrier Filter



PTFE Technology

