




- Made entirely from AISI 304 stainless steel for maximum durability and corrosion resistance
- Available in multiple sizes, designed to fit perfectly into our sterilization containers
- Can be customized with silicone holders, stainless steel partitions, and laser-marked ID labels
- Additionally, we can also create customized baskets, ensuring the perfect fit for any procedure
- Design protected
- Handles lift when needed, retract automatically, or foldable at 90°
- Larger bottom and lid perforations for better cleaning
- Finer perforations on the sides for instrument protection
- Optional feet in various heights for improved drainage and safe stackability
- Perforations extend to corners for full steam access and efficient drainage



FDA CLEARED 

You can find the complete product catalog under <https://www.ermis-medizintechnik.de/download-1/>

