



- Integrated safety lid for secure storage
- 100% Aluminum & Stainless steel construction, no plastic
- Noise-dampened and heat-resistant silicone-insulated handles
- Stackable design with clear airflow
- Compatible with paper (single-use), PTFE (up to 1200 sterilization cycles) & textile filters (multi-cycle)
- Seven color-coded lids and handles for efficient categorization and clear identification
- Designed for optional placement of indicator or identification labels on both handle sides
- Structurally enhanced walls retain shape through repeated high-temperature exposure

