



- Made entirely from AISI 304 stainless steel for maximum durability and corrosion resistance
- Available in multiple sizes, designed to fit perfectly into our sterilization containers
- Can be customized with silicone holders, stainless steel partitions, and laser-marked ID labels
- Additionally, we can also create customized baskets, ensuring the perfect fit for any procedure
- Integrated handle lifts easily and folds smoothly; slide-open tabs allow quick manual access
- Larger bottom and lid perforations for better cleaning
- Finer perforations on the sides for instrument protection
- Optional feet in various heights for improved drainage and safe stackability
- A corner opening ensures water drains from the lowest point, speeding up the drying process
- Designed to support heavier instruments, offering increased stability due to the basket's robust construction

