



Perforated-Side-Plate-Baskets




- Made entirely from AISI 304 stainless steel for maximum durability and corrosion resistance
- Available in multiple sizes, designed to fit perfectly into our sterilization containers
- Can be customized with silicone holders, stainless steel partitions, and laser-marked ID labels
- Additionally, we can also create customized baskets, ensuring the perfect fit for any procedure
- Handles lift when needed, retract automatically, or foldable at 90°
- Hybrid-basket with cleaning-efficient wire base and lid
- Sturdy perforated side plates safeguard instruments
- Rotating closing system ensures secure locking and foldable to avoid disturbance
- Bottom and lid available in three wire types adaptable to cleaning and protection needs


Type A



Type B



Type C



FDA CLEARED 

You can find the complete product catalog under <https://www.ermis-medizintechnik.de/download-1/>

